

HSP90 / Heat Shock Protein 90 Rabbit anti-Plasmodium falciparum Polyclonal (C-Terminus) (PerCP) Antibody - LS-C229695 - LSBio	
CatalogID:	LS-C229695
Target:	HSP90 / Heat Shock Protein 90
Host	Rabbit
Clonality:	Polyclonal
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	Plasmodium falciparum
Antigen Type:	Recombinant protein
Immunogen:	Recombinant partial PfHsp90 C-terminal 193 amino acids.
Specificity:	This protein A Affinity purified antibody detects an 86kD protein, corresponding to P. falciparum Hsp90 on SDS PAGE immunoblot; in samples from P. falciparum origin.
Epitope:	C-Terminus
Reactivity:	Plasmodium falciparum
Non-Reactivity:	Human, Yeast, Dictyostelium
Purification:	Protein A purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Immunofluorescence, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	