

AHSA2 Rabbit anti-Human Polyclonal (PerCP) Antibody - LS-C229711 - LSBio

CatalogID:	LS-C229711
Target:	AHA1, activator of heat shock 90kDa protein ATPase homolog 2 (yeast) (AHSA2)
Synonyms:	AHSA2 Antibody, Hch1 Antibody
Host	AHSA2 antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	AHSA2 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	AHSA2 antibody was raised against full length Aha2.
Specificity:	Cross reacts with Aha1.
Reactivity:	Human, Mouse, Rat
Purification:	Protein A purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences