

SLC9A3 / NHE3 Rabbit anti-Rat Polyclonal (PerCP) Antibody - LS-C229728 - LSBio

| | |
|-----------------------------|--|
| CatalogID: | LS-C229728 |
| Target: | solute carrier family 9, subfamily A (NHE3, cation proton antiporter 3), member 3 (SLC9A3) |
| Synonyms: | SLC9A3 Antibody, Na ⁺ /H ⁺ exchanger 3 Antibody, NHE3 Antibody, Sodium/hydrogen exchanger 3 Antibody, Na ⁽⁺⁾ /H ⁽⁺⁾ exchanger 3 Antibody, NHE-3 Antibody |
| Family / Subfamily: | Transporter / Monovalent cation:proton antiporter-1 |
| Host | SLC9A3 antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Conjugations: | Peridinin-chlorophyll-protein (PerCP) |
| Immunogen Species: | SLC9A3 / NHE3 antibody was raised against Rat |
| Antigen Type: | Synthetic peptide |
| Immunogen: | SLC9A3 / NHE3 antibody was raised against synthetic peptide mapping to AA 809 to 831 of rat sequence. |
| Specificity: | Detects ~84kD. |
| Reactivity: | Rat, Mouse |
| Purification: | Affinity purified |
| Presentation: | PBS, 0.09% sodium azide, 50% glycerol. |
| Recommended Storage: | Store at -20°C. |
| Usage Summary: | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested. |
| Uses: | ICC, Western blot (Optimal dilution to be determined by the researcher) |
| Size: | 100 µg |
| Requested From: | Japan |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences