

SCNN1G / ENaC Gamma Rabbit anti-Rat Polyclonal (aa629-650) (PerCP) Antibody - LS-C229780 - LSBio	
CatalogID:	LS-C229780
Target:	sodium channel, non-voltage-gated 1, gamma subunit (SCNN1G)
Synonyms:	SCNN1G Antibody, BESC3 Antibody, ENaC gamma Antibody, ENaCgamma Antibody, Gamma-NaCH Antibody, Gamma hENaC Antibody, HENaC Antibody, PHA1 Antibody, RENaC Antibody, SCNEG Antibody, ENaC gamma subunit Antibody, ENaCg Antibody, Gamma-ENaC Antibody
Family / Subfamily:	Ion Channel / Sodium channel - amiloride-sensitive
Host	SCNN1G antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	SCNN1G / ENaC Gamma antibody was raised against Rat
Antigen Type:	Synthetic peptide
Immunogen:	SCNN1G / ENaC Gamma antibody was raised against produced against the C-terminal tail of rat gamma ENaC AA629-650 (antibody designation L550).
Specificity:	Detects ~83kD.
Epitope:	aa629-650
Reactivity:	Rat, Mouse
Purification:	Affinity purified
Presentation:	PBS, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Immunofluorescence, Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	