

AQP3 / Aquaporin 3 Rabbit anti-Rat Polyclonal (C-Terminus) (PerCP) Antibody - LS-C229813 - LSBio

<b>CatalogID:</b>	LS-C229813
<b>Target:</b>	aquaporin 3 (Gill blood group) (AQP3)
<b>Synonyms:</b>	AQP3 Antibody, Aquaporin 3 (GIL blood group) Antibody, Aquaglyceroporin-3 Antibody, AQP-3 Antibody, Aquaporin 3 Antibody, Aquaporin 3 (Gill blood group) Antibody, Aquaporin-3 Antibody, GIL Antibody
<b>Family / Subfamily:</b>	Transporter / Aquaporin
<b>Host</b>	AQP3 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Conjugations:</b>	Peridinin-chlorophyll-protein (PerCP)
<b>Immunogen Species:</b>	AQP3 / Aquaporin 3 antibody was raised against Rat
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	AQP3 / Aquaporin 3 antibody was raised against sequence N-CHLEQPPPSTEAEENVKLAHMKHKEQI (Rat). Percent identity by BLAST analysis: Rat (100%).
<b>Specificity:</b>	Detects ~31.5kD. May detect larger glycosylated bands ~35-50kD.
<b>Epitope:</b>	C-Terminus
<b>Reactivity:</b>	Rat, Human, Mouse
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Immunofluorescence, Western blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences