

Amyloid Fibrils Rabbit anti-Human Polyclonal (PerCP) Antibody - LS-C229844 - LSBio

CatalogID:	LS-C229844
Target:	Amyloid Fibrils
Host	Amyloid Fibrils antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	Amyloid Fibrils antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	Amyloid Fibrils antibody was raised against fibrils prepared from human A beta 42 peptide.
Specificity:	Recognizes generic epitopes common to many amyloid fibrils and fibrillar oligomers, but not prefibrillar oligomers or native folded proteins.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Immunoprecipitation, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences