

mGluR1 + mGluR5 Mouse anti-Rat Monoclonal (aa824-1203) (PerCP) (S75-3) Antibody - LS-C230026 - LSBio	
CatalogID:	LS-C230026
Target:	Metabotropic Glutamate Receptor 1+5 (GRM1+GRM5)
Host	Metabotropic Glutamate Receptor 1+5 (GRM1+GRM5) antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	S75-3
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	mGluR1 + mGluR5 antibody was raised against Rat
Antigen Type:	Fusion protein
Immunogen:	mGluR1 + mGluR5 antibody was raised against fusion protein amino acids 824-1203 (cytoplasmic Cterminus) of rat mGluR5b, accession number P31424.
Specificity:	Detects ~130kD (both mGluR1 and mGluR5).
Epitope:	aa824-1203
Reactivity:	Rat, Human, Mouse
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	