

ACCN2 / ASIC Mouse anti-Mouse Monoclonal (aa460-526) (PerCP) (S271-44) Antibody - LS-C230061 - LSBio

CatalogID:	LS-C230061
Target:	acid-sensing (proton-gated) ion channel 1 (ASIC1)
Synonyms:	ASIC1 Antibody, Acid-sensing ion channel 1 Antibody, ACCN2 Antibody, Acid-sensing ion channel Antibody, ASIC1A Antibody, BNaC2 Antibody, Brain sodium channel 2 Antibody, HBNaC2 Antibody, ASIC Antibody, ASIC1B Antibody
Family / Subfamily:	Ion Channel / Sodium channel - amiloride-sensitive
Host	ASIC1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	S271-44
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	ACCN2 / ASIC antibody was raised against Mouse
Antigen Type:	Fusion protein
Immunogen:	ACCN2 / ASIC antibody was raised against fusion protein amino acids 460-526 (Cytoplasmic Cterminus) of mouse ASIC1, accession number Q6NXK8.
Specificity:	Detects ~60kD.
Epitope:	aa460-526
Reactivity:	Mouse, Rat
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences