

HSPH1 / HSP110 Rabbit anti-Hamster Polyclonal (PerCP) Antibody - LS-C230268 - LSBio	
CatalogID:	LS-C230268
Target:	heat shock 105kDa/110kDa protein 1 (HSPH1)
Synonyms:	HSPH1 Antibody, Heat shock 105kD beta Antibody, Heat shock 110 kDa protein Antibody, Heat shock protein 105 kDa Antibody, HSP105 Antibody, HSP110 Antibody, Heat-shock protein 105 kd Antibody, HSP105A Antibody, HSP105B Antibody, NY-CO-25 Antibody, KIAA0201 Antibody, Antigen NY-CO-25 Antibody, Heat shock 105kD alpha Antibody, Heat shock 105kDa protein 1 Antibody, Heat shock protein 110 Antibody
Host	HSPH1 antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	HSPH1 / HSP110 antibody was raised against Hamster
Antigen Type:	Synthetic peptide
Immunogen:	HSPH1 / HSP110 antibody was raised against synthetic peptide derived from the sequence of hamster Hsp110; sequence identical to human and mouse.
Specificity:	Detects ~110kD.
Reactivity:	Hamster, Human, Monkey, Mouse, Rat, Bovine, Sheep, Yeast
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014 © 2014 LifeSpan BioSciences	