

TLR4 Rabbit anti-Human Polyclonal (aa420-435) (RPE) Antibody - LS-C230289 - LSBio

CatalogID:	LS-C230289
Target:	toll-like receptor 4 (TLR4)
Synonyms:	TLR4 Antibody, CD284 Antibody, CD284 antigen Antibody, Homolog of Drosophila toll Antibody, TLR-4 Antibody, TOLL Antibody, ARMD10 Antibody, Toll-like receptor 4 Antibody, HToll Antibody
Family / Subfamily:	Toll-like Receptor / not assigned-Toll-like Receptor
Host	TLR4 antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	TLR4 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	TLR4 antibody was raised against developed against a synthetic peptide corresponding to amino acids 420-435 of human TLR4.
Specificity:	Detects ~75-80kD when tested against partial recombinant mouse TLR4 (extra-cellular portion plus His-tag).
Epitope:	aa420-435
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, 0.05% BSA, 0.05% sodium azide.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences