

SCNN1A / ENaC Apha Rabbit anti-Rat Polyclonal (aa46-68) (PerCP) Antibody - LS-C230320 - LSBio

<b>CatalogID:</b>	LS-C230320
<b>Target:</b>	sodium channel, non-voltage-gated 1 alpha subunit (SCNN1A)
<b>Synonyms:</b>	SCNN1A Antibody, Alpha hENaC Antibody, Alpha-NaCH Antibody, Alpha-ENaC Antibody, ENaCalpha Antibody, ENaCa Antibody, ENaC alpha Antibody, HENaC alpha Antibody, Alpha ENaC-2 Antibody, SCNEA Antibody, SCNN1 Antibody, BESC2 Antibody
<b>Family / Subfamily:</b>	Ion Channel / Sodium channel - amiloride-sensitive
<b>Host</b>	SCNN1A antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Conjugations:</b>	Peridinin-chlorophyll-protein (PerCP)
<b>Immunogen Species:</b>	SCNN1A / ENaC Apha antibody was raised against Rat
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	SCNN1A / ENaC Apha antibody was raised against produced against a synthetic peptide mapping to the N-terminal of the alpha subunit (amino acids 46-68) of rat Alpha ENaC (antibody designation 3560-2).
<b>Specificity:</b>	Detects ~85kD.
<b>Epitope:</b>	aa46-68
<b>Reactivity:</b>	Rat, Mouse
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, 0.09% sodium azide, 50% glycerol.
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Immunofluorescence, Western blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences