

TXNDC5 / ERP46 Mouse anti-Human Monoclonal (aa52-63) (PerCP) (2E7/7) Antibody - LS-C230362 - LSBio

CatalogID:	LS-C230362
Target:	thioredoxin domain containing 5 (endoplasmic reticulum) (TXNDC5)
Synonyms:	TXNDC5 Antibody, ENDOPDI Antibody, ER protein 46 Antibody, ERP46 Antibody, PDIA15 Antibody, STRF8 Antibody, TLP46 Antibody, UNQ364 Antibody, HCC-2 Antibody, Thioredoxin related protein Antibody, Thioredoxin-like protein p46 Antibody
Host	TXNDC5 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	2E7/7
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	TXNDC5 / ERP46 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	TXNDC5 / ERP46 antibody was raised against synthesized peptide - 12 amino acids long (ADGEDGQDPHSK) corresponding to residues 52 - 63 of the EndoPDI protein. Percent identity by BLAST analysis: Human, Chimpanzee (100%).
Specificity:	Detects a protein with a mass of ~48kD corresponding to EndoPDI.
Epitope:	aa52-63
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Immunofluorescence, Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences