

HSP90B1 / GP96 / GRP94 Rabbit anti-Mouse Polyclonal (C-Terminus) (FITC) Antibody - LS-C230393 - LSBio

CatalogID:	LS-C230393
Target:	heat shock protein 90kDa beta (Grp94), member 1 (HSP90B1)
Synonyms:	HSP90B1 Antibody, ECGP Antibody, Endoplasmin Antibody, gp96 homolog Antibody, gp96 Antibody, GRP-94 Antibody, GRP94 Antibody, Tumor rejection antigen 1 Antibody, TRA1 Antibody
Host	HSP90B1 antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	Fluorescein (FITC)
Immunogen Species:	HSP90B1 / GP96 / GRP94 antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	HSP90B1 / GP96 / GRP94 antibody was raised against synthetic peptide corresponding to the sequence near the C-terminus of mouse GRP94.
Specificity:	Detects ~94kD.
Epitope:	C-Terminus
Reactivity:	Mouse, Human, Rat, Bovine
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences