

PPP3CA / CCN1 / Calcineurin A Rabbit anti-Human Polyclonal (aa364-283) (RPE) Antibody - LS-C230431 - LSBio	
CatalogID:	LS-C230431
Target:	protein phosphatase 3, catalytic subunit, alpha isozyme (PPP3CA)
Synonyms:	PPP3CA Antibody, Calcineurin A Antibody, CALN Antibody, CALNA Antibody, CALNA1 Antibody, CNA1 Antibody, PPP2B Antibody, Calcineurin A alpha Antibody, CAM-PRP catalytic subunit Antibody, CCN1 Antibody, CNA Antibody, CNA alpha Antibody, pp2b-alpha Antibody
Family / Subfamily:	Protein Phosphatase / PP2B
Host	PPP3CA antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	PPP3CA / CCN1 / Calcineurin A antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PPP3CA / CCN1 / Calcineurin A antibody was raised against human Calcineurin A peptide (AA 364-283).
Specificity:	Identifies a band at ~61kD on a WB.
Epitope:	aa364-283
Reactivity:	Human, Mouse, Rat, Dog, Rabbit
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.4; 50% glycerol, 0.09% sodium azide.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	