

| DLG4 / PSD95 Mouse anti-Rat Monoclonal (FITC) (6G6) Antibody - LS-C230465 - LSBio | |
|---|---|
| CatalogID: | LS-C230465 |
| Target: | discs, large homolog 4 (Drosophila) (DLG4) |
| Synonyms: | DLG4 Antibody, Disks large homolog 4 Antibody, Presynaptic density protein 95 Antibody, PSD95 Antibody, SAP-90 Antibody, Synapse-associated protein 90 Antibody, PSD-95 Antibody, SAP90 Antibody, Tax interaction protein 15 Antibody, Discs large homolog 4 Antibody |
| Family / Subfamily: | MAGUK |
| Host | DLG4 antibody was produced in Mouse |
| Clonality: | Monoclonal |
| Isotype: | IgG2a,k |
| Clone Name: | 6G6 |
| Conjugations: | Fluorescein (FITC) |
| Immunogen Species: | DLG4 / PSD95 antibody was raised against Rat |
| Antigen Type: | Recombinant protein |
| Immunogen: | DLG4 / PSD95 antibody was raised against recombinant rat PSD-95. |
| Specificity: | Detects a ~100kD protein corresponding to the molecular mass of PSD-95 on SDS PAGE immunoblots. An additional protein of >100kD is also detected. Additional cross-reactive bands are detected at ~75kD and 50kD in rat and mouse samples. |
| Reactivity: | Rat, Human, Mouse, Bovine |
| Purification: | Protein G purified |
| Presentation: | PBS, pH 7.4, 0.09% sodium azide, 50% glycerol. |
| Recommended Storage: | Store at -20°C. |
| Usage Summary: | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested. |
| Uses: | IHC, ICC, Immunofluorescence, Western blot (Optimal dilution to be determined by the researcher) |
| Size: | 100 µg |
| Requested From: | Japan |
| Laboratory Reagent For In Vitro Research Use Only | |
| Not for resale without prior written consent from LifeSpan BioSciences, Inc. | |
| Created on 9/28/2014 | |
| © 2014 LifeSpan BioSciences | |