

DLG4 / PSD95 Mouse anti-Rat Monoclonal (PerCP) (6G6) Antibody - LS-C230468 - LSBio

<b>CatalogID:</b>	LS-C230468
<b>Target:</b>	discs, large homolog 4 (Drosophila) (DLG4)
<b>Synonyms:</b>	DLG4 Antibody, Disks large homolog 4 Antibody, Presynaptic density protein 95 Antibody, PSD95 Antibody, SAP-90 Antibody, Synapse-associated protein 90 Antibody, PSD-95 Antibody, SAP90 Antibody, Tax interaction protein 15 Antibody, Discs large homolog 4 Antibody
<b>Family / Subfamily:</b>	MAGUK
<b>Host</b>	DLG4 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a,k
<b>Clone Name:</b>	6G6
<b>Conjugations:</b>	Peridinin-chlorophyll-protein (PerCP)
<b>Immunogen Species:</b>	DLG4 / PSD95 antibody was raised against Rat
<b>Antigen Type:</b>	Recombinant protein
<b>Immunogen:</b>	DLG4 / PSD95 antibody was raised against recombinant rat PSD-95.
<b>Specificity:</b>	Detects a ~100kD protein corresponding to the molecular mass of PSD-95 on SDS PAGE immunoblots. An additional protein of >100kD is also detected. Additional cross-reactive bands are detected at ~75kD and 50kD in rat and mouse samples.
<b>Reactivity:</b>	Rat, Human, Mouse, Bovine
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, ICC, Immunofluorescence, Western blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences