

DLG4 / PSD95 Mouse anti-Rat Monoclonal (RPE) (6G6) Antibody - LS-C230469 - LSBio	
CatalogID:	LS-C230469
Target:	discs, large homolog 4 (Drosophila) (DLG4)
Synonyms:	DLG4 Antibody, Disks large homolog 4 Antibody, Presynaptic density protein 95 Antibody, PSD95 Antibody, SAP-90 Antibody, Synapse-associated protein 90 Antibody, PSD-95 Antibody, SAP90 Antibody, Tax interaction protein 15 Antibody, Discs large homolog 4 Antibody
Family / Subfamily:	MAGUK
Host	DLG4 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	6G6
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	DLG4 / PSD95 antibody was raised against Rat
Antigen Type:	Recombinant protein
Immunogen:	DLG4 / PSD95 antibody was raised against recombinant rat PSD-95.
Specificity:	Detects a ~100kD protein corresponding to the molecular mass of PSD-95 on SDS PAGE immunoblots. An additional protein of >100kD is also detected. Additional cross-reactive bands are detected at ~75kD and 50kD in rat and mouse samples.
Reactivity:	Rat, Human, Mouse, Bovine
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Immunofluorescence, Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	