

DNAJB1 / Hsp40 Mouse anti-Human Monoclonal (PerCP) (3B9.E6) Antibody - LS-C230519 - LSBio	
CatalogID:	LS-C230519
Target:	DnaJ (Hsp40) homolog, subfamily B, member 1 (DNAJB1)
Synonyms:	DNAJB1 Antibody, DnaJ protein homolog 1 Antibody, DNAJ1 Antibody, HDj-1 Antibody, Hdj1 Antibody, Heat shock 40 kDa protein 1 Antibody, Hsp40 Antibody, HSPF1 Antibody, RSPH16B Antibody, Sis1 Antibody, Radial spoke 16 homolog B Antibody, Heat shock 40kD protein 1 Antibody, Heat shock protein 40 Antibody, Human DnaJ protein 1 Antibody
Host	DNAJB1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	3B9.E6
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	DNAJB1 / Hsp40 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	DNAJB1 / Hsp40 antibody was raised against recombinant protein HSP40 (HDJ1).
Specificity:	Detects ~40kD. Does not cross-react with Hdj2 or YDJ1.
Reactivity:	Human, Mouse, Rat
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Immunofluorescence, Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014 © 2014 LifeSpan BioSciences	