

HSP90 / Heat Shock Protein 90 Mouse anti-Human Monoclonal (PerCP) (4F3.E8) Antibody - LS-C230536 - LSBio

CatalogID:	LS-C230536
Target:	HSP90 / Heat Shock Protein 90
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	4F3.E8
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	Human
Antigen Type:	Recombinant protein
Immunogen:	Recombinant Human Hsp90 purified from E. coli.
Specificity:	This antibody detects both alpha and beta forms of Hsp90 equally well.
Reactivity:	Human, Mouse, Rat, Wheat
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Immunofluorescence, Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences