

PTGES3 Mouse anti-Human Monoclonal (RPE) (JJ6) Antibody - LS-C230554 - LSBio

CatalogID:	LS-C230554
Target:	prostaglandin E synthase 3 (cytosolic) (PTGES3)
Synonyms:	PTGES3 Antibody, CPGES Antibody, p23 Antibody, Prostaglandin E synthase 3 Antibody, Telomerase-binding protein p23 Antibody, Hsp90 co-chaperone Antibody, TEBP Antibody
Host	PTGES3 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	JJ6
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	PTGES3 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	PTGES3 antibody was raised against recombinant Human full length p23 protein.
Specificity:	Detects ~23kD.
Reactivity:	Human, Mouse, Guinea pig, Rabbit, Chicken, Yeast
Purification:	Protein G purified
Presentation:	PBS, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences