

PDIA3 / ERp57 Mouse anti-Human Monoclonal (RPE) (Map.ERP57) Antibody - LS-C230586 - LSBio

CatalogID:	LS-C230586
Target:	protein disulfide isomerase family A, member 3 (PDIA3)
Synonyms:	PDIA3 Antibody, 58 kDa microsomal protein Antibody, ER protein 57 Antibody, ERp60 Antibody, ERp61 Antibody, Disulfide isomerase ER-60 Antibody, Endoplasmic reticulum P58 Antibody, GRP58 Antibody, HsT17083 Antibody, p58 Antibody, Phospholipase C-alpha Antibody, PI-PLC Antibody, ER protein 60 Antibody, ER60 Antibody, ERp57 Antibody, GRP57 Antibody, Protein disulfide-isomerase A3 Antibody
Host	PDIA3 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	Map.ERP57
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	PDIA3 / ERp57 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	PDIA3 / ERp57 antibody was raised against human Recombinant Erp57 (Grp58).
Specificity:	Detects ~57kD.
Reactivity:	Human, Monkey, Mouse, Rat, Bovine, Dog, Guinea pig, Hamster, Pig, Rabbit
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	200 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences