

HCN2 Mouse anti-Rat Monoclonal (aa761-863) (PerCP) (S71-37) Antibody - LS-C230654 - LSBio	
CatalogID:	LS-C230654
Target:	hyperpolarization activated cyclic nucleotide-gated potassium channel 2 (HCN2)
Synonyms:	HCN2 Antibody, BCNG2 Antibody, BCNG-2 Antibody, HAC-1 Antibody, HAC1 Antibody
Family / Subfamily:	Ion Channel / Potassium channel - cyclic nucleotide gated
Host	HCN2 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	S71-37
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	HCN2 antibody was raised against Rat
Antigen Type:	Fusion protein
Immunogen:	HCN2 antibody was raised against fusion protein amino acids 761-863 of rat HCN2, accession number Q9JKA9.
Specificity:	Detects ~95kD. No cross-reactivity against HCN1.
Epitope:	aa761-863
Reactivity:	Rat, Mouse
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Immunofluorescence, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	