

HSPD1 / HSP60 Rabbit anti-Plasmodium falciparum Polyclonal (PerCP) Antibody - LS-C230987 - LSBio

CatalogID:	LS-C230987
Target:	heat shock 60kDa protein 1 (chaperonin) (HSPD1)
Synonyms:	HSPD1 Antibody, 60 kDa chaperonin Antibody, CPN60 Antibody, Heat shock protein 65 Antibody, Heat shock protein 60 Antibody, HSP-60 Antibody, HSP60 Antibody, HuCHA60 Antibody, HLD4 Antibody, HSP65 Antibody, p60 lymphocyte protein Antibody, SPG13 Antibody, Chaperonin 60 Antibody, GROEL Antibody
Host	HSPD1 antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	HSPD1 / HSP60 antibody was raised against Plasmodium falciparum
Antigen Type:	Recombinant protein
Immunogen:	HSPD1 / HSP60 antibody was raised against recombinant full length PfHsp60.
Specificity:	This protein A Affinity purified antibody detects an 60kD protein, corresponding to P. falciparum Hsp60 on SDS PAGE immunoblot; in samples from P. falciparum origin.
Reactivity:	Plasmodium falciparum
Purification:	Protein A purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Immunofluorescence, Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences