

HSP70 / Heat Shock Protein 70 Mouse anti-Human Monoclonal (PerCP) (N27F3-4) Antibody - LS-C231029 - LSBio	
<b>CatalogID:</b>	LS-C231029
<b>Target:</b>	HSP70 / Heat Shock Protein 70
<b>Host</b>	Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	N27F3-4
<b>Conjugations:</b>	Peridinin-chlorophyll-protein (PerCP)
<b>Immunogen Species:</b>	Human
<b>Antigen Type:</b>	Recombinant protein
<b>Immunogen:</b>	Recombinant hsp70/hsc70.
<b>Specificity:</b>	Detects 72 and 73kD proteins corresponding to the molecular mass of inducible Hsp and Hsc on SDS PAGE immunoblots.
<b>Reactivity:</b>	Human, Monkey, Mouse, Rat, Bovine, Dog, Guinea pig, Hamster, Pig, Rabbit, Sheep, Chicken, Xenopus, Fish, Drosophila, Beluga, C. elegans, Pea
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide, 50% glycerol.
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, Immunoprecipitation, Flow Cytometry, Electron Microscopy (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µg
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	