

HSPD1 / HSP60 Mouse anti-Human Monoclonal (aa356-393) (RPE) (LK2) Antibody - LS-C231051 - LSBio	
<b>CatalogID:</b>	LS-C231051
<b>Target:</b>	heat shock 60kDa protein 1 (chaperonin) (HSPD1)
<b>Synonyms:</b>	HSPD1 Antibody, 60 kDa chaperonin Antibody, CPN60 Antibody, Heat shock protein 65 Antibody, Heat shock protein 60 Antibody, HSP-60 Antibody, HSP60 Antibody, HuCHA60 Antibody, HLD4 Antibody, HSP65 Antibody, p60 lymphocyte protein Antibody, SPG13 Antibody, Chaperonin 60 Antibody, GROEL Antibody
<b>Host</b>	HSPD1 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	LK2
<b>Conjugations:</b>	R. Phycoerythrin (RPE)
<b>Immunogen Species:</b>	HSPD1 / HSP60 antibody was raised against Human
<b>Antigen Type:</b>	Recombinant protein
<b>Immunogen:</b>	HSPD1 / HSP60 antibody was raised against recombinant human Hsp60, epitope aa356-393 of human protein.
<b>Specificity:</b>	Detects an ~60kD protein corresponding to the molecular mass of Hsp60 on SDS PAGE immunoblots.
<b>Epitope:</b>	aa356-393
<b>Reactivity:</b>	Human, Monkey, Mouse, Rat, Bovine, Dog, Guinea pig, Hamster, Pig, Rabbit, Chicken, Yeast, Trichina worm, E. coli, Helicobacter pylori, Salmonella typhimurium, Spinach, Whitefly
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, 0.09% sodium azide, 50% glycerol.
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µg
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	