

HSF2 Rat anti-Mouse Monoclonal (PerCP) (3.00E+02) Antibody - LS-C231062 - LSBio

CatalogID:	LS-C231062
Target:	heat shock transcription factor 2 (HSF2)
Synonyms:	HSF2 Antibody, HSTF 2 Antibody, Heat shock factor protein 2 Antibody, HSF 2 Antibody, HSTF2 Antibody
Host	HSF2 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	3.00E+02
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	HSF2 antibody was raised against Mouse
Antigen Type:	Recombinant protein
Immunogen:	HSF2 antibody was raised against purified recombinant mouse HSF2 protein.
Specificity:	Detects an ~69kD protein in unstressed cells corresponding to the molecular mass HSF2 on SDS PAGE immunoblots.
Reactivity:	Mouse, Human, Monkey, Rat, Bovine, Dog, Guinea pig, Hamster, Pig, Rabbit, Sheep
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences