

Nitrotyrosine Mouse Monoclonal (PerCP) (39B6) Antibody - LS-C231161 - LSBio

CatalogID:	LS-C231161
Target:	Nitrotyrosine
Host	Nitrotyrosine antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	39B6
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Antigen Type:	Hapten
Immunogen:	Nitrotyrosine antibody was raised against nitrotyrosine.
Specificity:	Recognizes 3-nitrotyrosine moieties. No detectable cross-reactivity with non-nitrated tyrosine. Not species specific.
Purification:	Protein G purified
Presentation:	PBS, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Immunofluorescence, Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences