

HSP90 / Heat Shock Protein 90 Mouse anti-Human Monoclonal (RPE) (8D3) Antibody - LS-C231195 - LSBio

CatalogID:	LS-C231195
Target:	HSP90 / Heat Shock Protein 90
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgM
Clone Name:	8D3
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	Human
Antigen Type:	Purified protein
Immunogen:	Ah receptor (Aryl hydrocarbon receptor).
Specificity:	Immunoprecipitates 90kD proteins corresponding to the molecular mass of Hsp90. Co-immunoprecipitates Hsp90 complexes, including Hsp70, Hop, Ah receptors, glucocorticoid receptors, heme-regulated eukaryotic initiation factor 2 alpha (eIF-2 alpha) kinase (HRI).
Reactivity:	Human, Mouse, Rat, Rabbit
Purification:	PEG purified
Presentation:	PBS, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	200 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences