

HSPB2 / HSP27 Mouse anti-Human Monoclonal (PerCP) (8A7) Antibody - LS-C231212 - LSBio	
CatalogID:	LS-C231212
Target:	heat shock 27kDa protein 2 (HSPB2)
Synonyms:	HSPB2 Antibody, DMPK-binding protein Antibody, Heat shock protein beta-2 Antibody, HSP27 Antibody, Heat-shock protein beta-2 Antibody, Hs.78846 Antibody, MKBP Antibody, Heat shock 27kD protein 2 Antibody, Heat shock 27kDa protein 2 Antibody, LOH11CR1K Antibody
Host	HSPB2 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG
Clone Name:	8A7
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	HSPB2 / HSP27 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	HSPB2 / HSP27 antibody was raised against hsp27 peptide.
Specificity:	Detects an 25kD or 27kD protein corresponding to the molecular mass of Hsp25/27 on SDS PAGE immunoblots.
Reactivity:	Human, Mouse, Rat, Bovine, Dog, Guinea pig, Hamster
Purification:	Protein G purified
Presentation:	PBS, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Immunofluorescence, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	