

CRYAA / Alpha A Crystallin Mouse anti-Human Monoclonal (PerCP) (1H3.B8) Antibody - LS-C231246 - LSBio

CatalogID:	LS-C231246
Target:	crystallin, alpha A (CRYAA)
Synonyms:	CRYAA Antibody, Alpha-crystallin A chain Antibody, Crystallin, alpha A Antibody, CRYA1 Antibody, Crystallin, alpha-1 Antibody, Heat shock protein beta-4 Antibody, HSPB4 Antibody
Host	CRYAA antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	1H3.B8
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	CRYAA / Alpha A Crystallin antibody was raised against Human
Antigen Type:	Native protein
Immunogen:	CRYAA / Alpha A Crystallin antibody was raised against native Alpha Crystallin.
Specificity:	Detects alpha -crystallin at ~20kD. Does not cross-react with alpha B-crystallin, beta L-crystallin, [Character 392]H- crystallin, gamma -crystallin, Hsp25, Hsp27 or Hsp47 proteins.
Reactivity:	Human, Mouse, Rat, Bovine
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, in 50% glycerol, 0.09% sodium azide
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences