

HSP90AA / HSP90 Alpha Mouse anti-Human Monoclonal (FITC) (2G5.G3) Antibody - LS-C231260 - LSBio

CatalogID:	LS-C231260
Target:	HSP90 Alpha (HSP90AA)
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Clone Name:	2G5.G3
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Human
Antigen Type:	Purified protein
Immunogen:	Human Hsp90alpha.
Specificity:	Hsp90 alpha -specific (>96% alpha -specific by ELISA).
Reactivity:	Human, Mouse, Rat
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Immunofluorescence, Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences