

HSP90AA / HSP90 Alpha Mouse anti-Human Monoclonal (PerCP) (2G5.G3) Antibody - LS-C231263 - LSBio

<b>CatalogID:</b>	LS-C231263
<b>Target:</b>	HSP90 Alpha (HSP90AA)
<b>Host</b>	Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1,k
<b>Clone Name:</b>	2G5.G3
<b>Conjugations:</b>	Peridinin-chlorophyll-protein (PerCP)
<b>Immunogen Species:</b>	Human
<b>Antigen Type:</b>	Purified protein
<b>Immunogen:</b>	Human Hsp90alpha.
<b>Specificity:</b>	Hsp90 alpha -specific (>96% alpha -specific by ELISA).
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide, 50% glycerol.
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, ICC, Immunofluorescence, Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences