

HSPA8 / HSC70 Mouse anti-Human Monoclonal (RPE) (1F2-H5) Antibody - LS-C231281 - LSBio	
CatalogID:	LS-C231281
Target:	heat shock 70kDa protein 8 (HSPA8)
Synonyms:	HSPA8 Antibody, Heat shock 70 kDa protein 8 Antibody, Heat shock 70kd protein 10 Antibody, Heat shock 70kD protein 8 Antibody, Heat shock cognate protein 54 Antibody, HSP71 Antibody, HSPA10 Antibody, HSC71 Antibody, HSP73 Antibody, LPS-associated protein 1 Antibody, Heat shock 70kDa protein 8 Antibody, NIP71 Antibody, HSC54 Antibody, HSC70 Antibody, LAP1 Antibody
Host	HSPA8 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	1F2-H5
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	HSPA8 / HSC70 antibody was raised against Human
Antigen Type:	Purified protein
Immunogen:	HSPA8 / HSC70 antibody was raised against human Hsc70.
Specificity:	Detects ~73kD. Does not cross react with Hsp70.
Reactivity:	Human, Mouse, Rat
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014 © 2014 LifeSpan BioSciences	