

Ubiquitin Mouse anti-Bovine Monoclonal (PerCP) (5B9-B3) Antibody - LS-C231314 - LSBio	
<b>CatalogID:</b>	LS-C231314
<b>Target:</b>	Ubiquitin
<b>Host</b>	Ubiquitin antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a,k
<b>Clone Name:</b>	5B9-B3
<b>Conjugations:</b>	Peridinin-chlorophyll-protein (PerCP)
<b>Immunogen Species:</b>	Ubiquitin antibody was raised against Bovine
<b>Antigen Type:</b>	Native protein
<b>Immunogen:</b>	Ubiquitin antibody was raised against native bovine ubiquitin, conjugated to KLH.
<b>Specificity:</b>	This antibody recognizes ~10kD corresponding to free ubiquitin.
<b>Reactivity:</b>	Bovine, Human, Mouse, Rat
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µg
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	