

Ubiquitin Mouse anti-Bovine Monoclonal (FITC) (6C11-B3) Antibody - LS-C231362 - LSBio

CatalogID:	LS-C231362
Target:	Ubiquitin
Host	Ubiquitin antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	6C11-B3
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Ubiquitin antibody was raised against Bovine
Antigen Type:	Native protein
Immunogen:	Ubiquitin antibody was raised against native bovine ubiquitin, conjugated to KLH.
Specificity:	This antibody recognizes ~10kD corresponding to free ubiquitin.
Reactivity:	Bovine, Human, Mouse, Rat
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences