

AHSA1 / AHA1 Rat anti-Mouse Monoclonal (PerCP) (25F2.D9) Antibody - LS-C231382 - LSBio	
CatalogID:	LS-C231382
Target:	AHA1, activator of heat shock 90kDa protein ATPase homolog 1 (yeast) (AHSA1)
Synonyms:	AHSA1 Antibody, AHA1 Antibody, C14orf3 Antibody, p38 Antibody
Host	AHSA1 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	25F2.D9
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	AHSA1 / AHA1 antibody was raised against Mouse
Antigen Type:	Purified protein
Immunogen:	AHSA1 / AHA1 antibody was raised against mouse Aha1.
Specificity:	Detects ~38kD. Can run up to 45kD on SDS Page.
Reactivity:	Mouse, Human, Rat
Purification:	Protein G purified
Presentation:	0.02 M potassium phosphate, 0.15 M sodium chloride, pH 7.2.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Immunofluorescence, Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	