

HIF1A / HIF1 Alpha Mouse anti-Mouse Monoclonal (aa329-530) (PerCP) (ESEE122) Antibody - LS-C231419 - LSBio	
CatalogID:	LS-C231419
Target:	hypoxia inducible factor 1, alpha subunit (basic helix-loop-helix transcription factor) (HIF1A)
Synonyms:	HIF1A Antibody, ARNT-interacting protein Antibody, ARNT interacting protein Antibody, BHLHe78 Antibody, HIF1 Alpha Antibody, Hypoxia-inducible factor1alpha Antibody, HIF-1A Antibody, HIF1 Antibody, MOP1 Antibody, Member of PAS superfamily 1 Antibody, PASD8 Antibody, HIF-1-alpha Antibody, HIF-1alpha Antibody, HIF1-ALPHA Antibody, Member of PAS protein 1 Antibody
Host	HIF1A antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	ESEE122
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	HIF1A / HIF1 Alpha antibody was raised against Mouse
Antigen Type:	Recombinant protein
Immunogen:	HIF1A / HIF1 Alpha antibody was raised against recombinant fragment corresponding to amino acids 329-530
Specificity:	Specific for HIF1Alpha.
Epitope:	aa329-530
Reactivity:	Mouse, Human, Rat, Bovine
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC - Paraffin, IHC - Frozen, ICC, Immunofluorescence, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	