

HSPB8 / H11 / HSP22 Mouse anti-Human Monoclonal (PerCP) (3C12-H11) Antibody - LS-C231432 - LSBio	
CatalogID:	LS-C231432
Target:	heat shock 22kDa protein 8 (HSPB8)
Synonyms:	HSPB8 Antibody, Alpha-crystallin C chain Antibody, CMT2L Antibody, Crystallin alpha c Antibody, DHMN2 Antibody, E2-induced gene 1 protein Antibody, CRYAC Antibody, E2IG1 Antibody, HMN2 Antibody, H11 Antibody, Heat shock 27kDa protein 8 Antibody, Heat shock protein beta-8 Antibody, Heat shock 22kDa protein 8 Antibody, HMN2A Antibody, HSP22 Antibody, Protein kinase H11 Antibody
Family / Subfamily:	Protein Kinase / H11
Host	HSPB8 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Clone Name:	3C12-H11
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	HSPB8 / H11 / HSP22 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	HSPB8 / H11 / HSP22 antibody was raised against his-tagged human recombinant hsp22
Specificity:	Detects endogenous and exogenous hsp22 in monomeric, dimeric and tetrameric forms in WB. Does not cross react with alpha crystallin. Exogenous hsp22 detected in dimeric form.
Reactivity:	Human, Mouse, Rat
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Immunofluorescence, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	