

MDC1 Mouse anti-Mouse Monoclonal (N-Terminus) (RPE) (P2B11) Antibody - LS-C231467 - LSBio	
<b>CatalogID:</b>	LS-C231467
<b>Target:</b>	mediator of DNA-damage checkpoint 1 (MDC1)
<b>Synonyms:</b>	MDC1 Antibody, Em:AB023051.5 Antibody, NFBD1 Antibody, KIAA0170 Antibody
<b>Host</b>	MDC1 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	P2B11
<b>Conjugations:</b>	R. Phycoerythrin (RPE)
<b>Immunogen Species:</b>	MDC1 antibody was raised against Mouse
<b>Antigen Type:</b>	Recombinant protein
<b>Immunogen:</b>	MDC1 antibody was raised against gST-tagged recombinant protein corresponding to mouse MDC1 at and around the N-terminus
<b>Specificity:</b>	Detects ~184kD. This antibody recognizes MDC1 at or around the N-terminus.
<b>Epitope:</b>	N-Terminus
<b>Reactivity:</b>	Mouse, Human, Chimpanzee, Bovine
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Immunofluorescence, Western blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	