

AHSA1 / AHA1 Rabbit anti-Mouse Polyclonal (PerCP) Antibody - LS-C231481 - LSBio

CatalogID:	LS-C231481
Target:	AHA1, activator of heat shock 90kDa protein ATPase homolog 1 (yeast) (AHSA1)
Synonyms:	AHSA1 Antibody, AHA1 Antibody, C14orf3 Antibody, p38 Antibody
Host	AHSA1 antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	AHSA1 / AHA1 antibody was raised against Mouse
Antigen Type:	Purified protein
Immunogen:	AHSA1 / AHA1 antibody was raised against mouse Aha1.
Specificity:	Detects ~38kD.
Reactivity:	Mouse, Human, Rat
Purification:	Protein A purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences