

Phosphotyrosine Mouse Monoclonal (PerCP) (G104) Antibody - LS-C231514 - LSBio

CatalogID:	LS-C231514
Target:	Phosphotyrosine
Host	Phosphotyrosine antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	G104
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Antigen Type:	Synthetic peptide
Immunogen:	Phosphotyrosine antibody was raised against phosphotyrosine, alanine and glycine in a 1:1:1 ratio polymerized in the presence of keyhole limpet hemocyanin with 1-ethyl-3-(3'-dimentrylamino)propyl carbodiimide.
Specificity:	Reacts with phosphotyrosine, and detects the presence of phosphotyrosine in both un-stimulated and stimulated cell lysates. Does not cross-react with phosphoserine or phosphothreonine.
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Immunofluorescence, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences