

DUX4 Mouse anti-Human Monoclonal (C-Terminus) (FITC) (P2B1) Antibody - LS-C231549 - LSBio	
CatalogID:	LS-C231549
Target:	double homeobox 4 (DUX4)
Synonyms:	DUX4 Antibody, Double homeobox protein 10 Antibody, Double homeobox protein DUX10 Antibody, Double homeobox 4 Antibody, Double homeobox protein 4/10 Antibody, Double homeobox protein 4 Antibody, DUX10 Antibody
Host	DUX4 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	P2B1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	DUX4 antibody was raised against Human
Antigen Type:	Fusion protein
Immunogen:	DUX4 antibody was raised against c-terminal 76 amino acids of DUX4 with glutathione-s-transferase (gst) tag
Specificity:	This antibody recognizes human DUX4. It does not cross-react with DUX4c.
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Immunofluorescence, Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014 © 2014 LifeSpan BioSciences	