

ABCC8 / SUR1 Mouse anti-Rat Monoclonal (aa1548-1582) (PerCP) (S289-16) Antibody - LS-C231739 - LSBio	
CatalogID:	LS-C231739
Target:	ATP-binding cassette, sub-family C (CFTR/MRP), member 8 (ABCC8)
Synonyms:	ABCC8 Antibody, ABC36 Antibody, HI Antibody, PHHI Antibody, Sulfonylurea receptor Antibody, SUR Antibody, Sulfonylurea receptor 1 Antibody, MRP8 Antibody, SUR1delta2 Antibody, HHF1 Antibody, HRINS Antibody, SUR1 Antibody, TNDM2 Antibody
Family / Subfamily:	Transporter / ATP-binding cassette - ABCC/MRP
Host	ABCC8 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	S289-16
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	ABCC8 / SUR1 antibody was raised against Rat
Antigen Type:	Fusion protein
Immunogen:	ABCC8 / SUR1 antibody was raised against fusion protein amino acids 1548-1582 (cytoplasmic C-terminus) of rat SUR1
Specificity:	Detects ~160kD. Does not cross-react with SUR2B.
Epitope:	aa1548-1582
Reactivity:	Rat, Mouse, Hamster
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	