

SOD1 / SOD Rabbit anti-Rat Polyclonal (PerCP) Antibody - LS-C231768 - LSBio

| | |
|-----------------------------|--|
| CatalogID: | LS-C231768 |
| Target: | superoxide dismutase 1, soluble (SOD1) |
| Synonyms: | SOD1 Antibody, ALS Antibody, ALS1 Antibody, Cu/Zn superoxide dismutase Antibody, HSod1 Antibody, Homodimer Antibody, IPOA Antibody, SOD, soluble Antibody, SOD Antibody, Superoxide dismutase [Cu-Zn] Antibody, Superoxide dismutase 1 Antibody, Indophenoloxidase A Antibody, Superoxide dismutase, cystolic Antibody |
| Host | SOD1 antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Conjugations: | Peridinin-chlorophyll-protein (PerCP) |
| Immunogen Species: | SOD1 / SOD antibody was raised against Rat |
| Antigen Type: | Purified protein |
| Immunogen: | SOD1 / SOD antibody was raised against rat Cu/Zn SOD. |
| Specificity: | Detects a ~23kD (human) and 19kD (other species) proteins corresponding to the molecular mass of Cu/Zn superoxide dismutase (SOD) on SDS PAGE immunoblots. |
| Reactivity: | Rat, Human, Mouse, Bovine |
| Purification: | Immunoaffinity purified |
| Presentation: | PBS, pH 7.0, 0.09% sodium azide, 50% glycerol |
| Recommended Storage: | Store at -20°C. |
| Usage Summary: | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested. |
| Uses: | IHC, Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher) |
| Size: | 100 µg |
| Requested From: | Japan |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences