

Amyloid Oligomers Rabbit anti-Human Polyclonal (PerCP) Antibody - LS-C231801 - LSBio

CatalogID:	LS-C231801
Target:	Amyloid Oligomers
Host	Rabbit
Clonality:	Polyclonal
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	Human
Antigen Type:	Synthetic peptide
Immunogen:	Synthetic molecular mimic of soluble oligomers.
Specificity:	Recognizes all types of amyloid oligomers. Appears to recognize a peptide backbone epitope that is common to amyloid oligomers, but is not found in native proteins, amyloidogenic monomer or mature amyloid fibrils.
Reactivity:	Human, Mouse, Rat, Eukaryote
Purification:	Protein A purified
Presentation:	PBS, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Immunofluorescence, Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences