

LAMP1 / CD107a Mouse anti-Rat Monoclonal (PerCP) (Ly1C6) Antibody - LS-C231834 - LSBio	
CatalogID:	LS-C231834
Target:	lysosomal-associated membrane protein 1 (LAMP1)
Synonyms:	LAMP1 Antibody, CD107a Antibody, CD107a antigen Antibody, LAMP-1 Antibody, LGP120 Antibody, LAMPA Antibody
Host	LAMP1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	Ly1C6
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	LAMP1 / CD107a antibody was raised against Rat
Antigen Type:	Membrane preparation
Immunogen:	LAMP1 / CD107a antibody was raised against rat liver lysosomal membrane preparations.
Specificity:	Detects an ~120kD protein corresponding to the molecular mass of LAMP1 on SDS PAGE immunoblots.
Reactivity:	Rat, Hamster
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ICC, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014 © 2014 LifeSpan BioSciences	