

HSPB1 / HSP27 Mouse anti-Human Monoclonal (PerCP) (5D12-A12) Antibody - LS-C231851 - LSBio

CatalogID:	LS-C231851
Target:	heat shock 27kDa protein 1 (HSPB1)
Synonyms:	HSPB1 Antibody, 28 kDa heat shock protein Antibody, Heat shock 27 kDa protein Antibody, Heat shock 27kD protein 1 Antibody, Heat shock 27kDa protein 1 Antibody, HSP27 Antibody, HSP28 Antibody, HSP 27 Antibody, Hsp25 Antibody, SRP27 Antibody, Stress-responsive protein 27 Antibody, CMT2F Antibody, Heat shock 27 kd protein Antibody, Heat shock protein beta-1 Antibody, HMN2B Antibody, HS.76067 Antibody
Host	HSPB1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b,k
Clone Name:	5D12-A12
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	HSPB1 / HSP27 antibody was raised against Human
Antigen Type:	Purified protein
Immunogen:	HSPB1 / HSP27 antibody was raised against human Hsp27.
Specificity:	Detects ~27kD proteins corresponding to Hsp27 on SDS Page Western Blots. Has no cross-reactivity to alphaB crystallin.
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences