

HSP90 / Heat Shock Protein 90 Mouse anti-Achlya ambisexualis Monoclonal (PerCP) (AC-16)
Antibody - LS-C231890 - LSBio

CatalogID:	LS-C231890
Target:	HSP90 / Heat Shock Protein 90
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	AC-16
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	Achlya ambisexualis
Antigen Type:	Purified protein
Immunogen:	Heat shock protein 900 from the water mold Achlya ambisexualis.
Specificity:	This antibody is reactive with both the constitutive and the inducible form of Hsp90. It does not bind to the native form and does not recognize Hsp90 from E. coli or yeast.
Reactivity:	Achlya ambisexualis, Human, Mouse, Rat, Rabbit, Chicken, Spodoptera frugiperda, Wheat Germ
Non-Reactivity:	Yeast, E. coli
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot (Optimal dilution to be determined by the researcher)
Size:	200 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences