

KDEL / KDEL Receptor Mouse anti-Bovine Monoclonal (PerCP) (KR-10) Antibody - LS-C231906 - LSBio

CatalogID:	LS-C231906
Target:	KDEL Receptor (KDEL)
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	KR-10
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	Bovine
Antigen Type:	Synthetic peptide - KLH conjugated
Immunogen:	A 21 residue synthetic peptide based on the bovine KDEL receptor and the peptide coupled to KLH.
Specificity:	Detects an ~25kD protein corresponding to the molecular mass of KDEL on SDS-PAGE immunoblots.
Reactivity:	Bovine, Human, Monkey, Mouse, Rat, Dog, Hamster, Pig, Rabbit, Sheep, Chicken, Xenopus, Drosophila
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Immunofluorescence, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences