

**HSP90 / Heat Shock Protein 90 Mouse anti-Chicken Monoclonal (PerCP) (D7A) Antibody - LS-C231923 - LSBio**

<b>CatalogID:</b>	LS-C231923
<b>Target:</b>	HSP90 / Heat Shock Protein 90
<b>Host</b>	Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	D7A
<b>Conjugations:</b>	Peridinin-chlorophyll-protein (PerCP)
<b>Immunogen Species:</b>	Chicken
<b>Antigen Type:</b>	Purified protein
<b>Immunogen:</b>	Full length protein Hsp90 purified from chicken brain.
<b>Specificity:</b>	Recognizes 90kD proteins corresponding to the molecular mass of Hsp90. Hsp90 alpha specific for human samples. Can isolate complexes of Hsp90, Src kinase and cec37.
<b>Reactivity:</b>	Chicken, Human, Mouse, Rat, Bovine, Pig, Rabbit
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, 0.09% sodium azide, 50% glycerol.
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences